

FAM84B Antibody



PACO09224

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FAM84B

Uniprot

Q96KN1

Synonyms:

family with sequence similarity 84, member B; FAM84B; BCMP101; NSE2 ;

Immunogen:

Human FAM84B.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A