

FAM96B Antibody



PACO09225

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FAM96B

Uniprot

Q9Y3D0

Synonyms:

family with sequence similarity 96, member B; FAM96B; CGI-128 ;

Immunogen:

Human FAM96B.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A