

FANCL Antibody



PACO09231

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FANCL

Uniprot

Q9NW38

Synonyms:

Fanconi anemia, complementation group L; FANCL; FAAP43; FLJ10335; PHF9; POG ;

Immunogen:

Human FANCL

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A