

FARP2 Antibody



PACO09233

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FARP2

Uniprot

O94887

Synonyms:

FERM, RhoGEF and pleckstrin domain protein 2; FARP2; FIR; FRG; KIAA0793; PLEKHC3 ;

Immunogen:

Human FARP2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A