

# FBXL14 Antibody



PACO09247

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

FBXL14

**Uniprot**

Q8N1E6

**Synonyms:**

F-box and leucine-rich repeat protein 14; FBXL14; Fbl14; MGC40195 ;

**Immunogen:**

Human FBXL14.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A