

# FBXO6 Antibody



PACO09257

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FBXO6

**Uniprot**

Q9NRD1

**Synonyms:**

F-box protein 6; FBXO6; FBG2; FBS2; FBX6; Fbx6b ;

**Immunogen:**

Human FBXO6.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A