

FBXO8 Antibody



PACO09259

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

FBXO8

Uniprot

Q9NRD0

Synonyms:

F-box protein 8; FBXO8; DC10; FBS; FBX8 ;

Immunogen:

Human FBXO8.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A