

FBXW2 Antibody



PACO09262

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FBXW2

Uniprot

Q9UKT8

Synonyms:

F-box and WD repeat domain containing 2; FBXW2; FBW2; Fwd2; MGC117371; Md6 ;

Immunogen:

Human FBXW2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A