

FCRLB Antibody



PACO09269

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FCRLB

Uniprot

Q6BAA4

Synonyms:

Fc receptor-like B; FCRLB; FCRL2; FCRLM2; FCRLY; FLJ31052; FREB-2; FcRY; RP11-474116.6 ;

Immunogen:

Human FCRLB.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A