

FEM1A Antibody



PACO09273

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FEM1A

Uniprot

Q9BSK4

Synonyms:

fem-1 homolog a (C. elegans); FEM1A; DKFZp762M136; EPRAP ;

Immunogen:

Human FEM1A.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A