

FGD3 Antibody



PACO09281

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

FGD3

Uniprot

Q5JSP0

Synonyms:

FYVE, RhoGEF and PH domain containing 3; FGD3; FLJ00004; MGC117260; ZFYVE5 ;

Immunogen:

Human FGD3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A