

FGFR1OP Antibody



PACO09287

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FGFR1OP

Uniprot

O95684

Synonyms:

FGFR1 oncogene partner; FGFR1OP; FOP ;

Immunogen:

Human FGFR1OP.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A