

FHL3 Antibody



PACO09297

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

FHL3

Uniprot

Q13643

Synonyms:

four and a half LIM domains 3; FHL3; MGC19547; MGC23614; MGC8696; SLIM2 ;

Immunogen:

Human FHL3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A