

# FIS1 Antibody



PACO09300

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FIS1

**Uniprot**

Q9Y3D6

**Synonyms:**

fission 1 (mitochondrial outer membrane) homolog (S. cerevisiae); FIS1; CGI-135; TTC11 ;

**Immunogen:**

Human FIS1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A