

# FMR1NB Antibody



PACO09327

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

FMR1NB

**Uniprot**

Q8N0W7

**Synonyms:**

fragile X mental retardation 1 neighbor; FMR1NB; CT37; FLJ25736; NY-SAR-35 ;

**Immunogen:**

Human FMR1NB.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A