

FNBP1L Antibody



PACO09330

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FNBP1L

Uniprot

Q5T0N5

Synonyms:

formin binding protein 1-like; FNBP1L; C1orf39; TOCA1 ;

Immunogen:

Human FNBP1L.

Storage:

0.1M NaHCO₃, 0.1M Glycine, 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A