

FOXO6 Antibody



PACO09341

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FOXO6

Uniprot

A8MYZ6

Synonyms:

forkhead box protein O6;

Immunogen:

Human FOXO6.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A