

FSCN3 Antibody



PACO09351

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

FSCN3

Uniprot

Q9NQT6

Synonyms:

fascin homolog 3, actin-bundling protein, testicular (Strongylocentrotus purpuratus);

Immunogen:

Human FSCN3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A