

FXYD4 Antibody



PACO09370

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FXYD4

Uniprot

P59646

Synonyms:

FXYD domain containing ion transport regulator 4; FXYD4; CHIF ;

Immunogen:

Human FXYD4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A