

FZR1 Antibody



PACO09374

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

FZR1

Uniprot

Q9UM11

Synonyms:

fizzy/cell division cycle 20 related 1 (Drosophila); FZR1; CDC20C; CDH1; FZR; FZR2; HCDH; HCDH1; KIAA1242 ;

Immunogen:

Human FZR1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A