

# GABBR1 Antibody



PACO09380

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:****Protein Background:****Gene ID:**

GABBR1

**Uniprot**

Q9UBS5

**Synonyms:**

gamma-aminobutyric acid, (GABA) B receptor, 1; GABBR1; FLJ92613; GABAB(1e); GABABR1; GABBR1-3; GPRC3A; dJ271M21.1.1; dJ271M21.1.2; hGB1a ;

**Immunogen:**

Human GABBR1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A