

GALNT15 Antibody



PACO09405

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GALNT15

Uniprot

Q8N3T1

Synonyms:

UDP-N-acetyl-alpha-D-galactosamine: polypeptide N-acetylgalactosaminyltransferase-like 2; GALNTL2; DKFZp686H1113; GALNT13; GALNT15; GALNT7; PIH5 ;

Immunogen:

Human GALNTL2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A