

PACO09408

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat, Zebrafish

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GAMT

Uniprot

Q14353

Synonyms:

guanidinoacetate N-methyltransferase; GAMT; PIG2; TP53I2 ;

Immunogen:

Human GAMT.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A