

GCC1 Antibody



PACO09430

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GCC1

Uniprot

Q96CN9

Synonyms:

GRIP and coiled-coil domain containing 1; GCC1; FLJ22035; GCC1P; GCC88; MGC20706 ;

Immunogen:

Human GCC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A