

GCNT2 Antibody



PACO09436

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GCNT2

Uniprot

Q06430

Synonyms:

glucosaminyl (N-acetyl) transferase 2, I-branching enzyme (I blood group); GCNT2; CCAT; GCNT2C; GCNT5; IGNT; II; MGC163396; NACGT1; NAGCT1; ULG3; bA360O19.2; bA421M1.1 ;

Immunogen:

Human GCNT2.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A