

GCOM1 Antibody



PACO09437

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

GCOM1

Uniprot

H8Y6P7

Synonyms:

GRINL1A complex locus; GCOM1; FLJ30973; GRINL1A; Gcom2; Gup1; Gup2; MGC126694; MGC138353 ;

Immunogen:

Human GCOM1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A