

# GDPD5 Antibody



PACO09444

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

GDPD5

**Uniprot**

Q8WTR4

**Synonyms:**

glycerophosphodiester phosphodiesterase domain containing 5; GDPD5; GDE2; PP1665  
;

**Immunogen:**

Human GDPD5.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20°C; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A