

GFI1 Antibody



PACO09451

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GFI1

Uniprot

Q99684

Synonyms:

growth factor independent 1 transcription repressor; GFI1; FLJ94509; GFI-1; ZNF163 ;

Immunogen:

Human GFI1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A