

# GIMAP7 Antibody



PACO09468

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IF

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GIMAP7

**Uniprot**

Q8NHV1

**Synonyms:**

GTPase, IMAP family member 7; GIMAP7; IAN7; MGC27027; hIAN7 ;

**Immunogen:**

Human GIMAP7.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A