

GJA3 Antibody



PACO09474

Product Information

Size:

50ul

Reactivity:

Human

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GJA3

Uniprot

Q9Y6H8

Synonyms:

gap junction protein, alpha 3, 46kDa; GJA3; CX46; CZP3 ;

Immunogen:

Human GJA3.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A