

# GLYATL1 Antibody



PACO09495

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GLYATL1

**Uniprot**

Q969I3

**Synonyms:**

glycine-N-acyltransferase-like 1; GLYATL1; FLJ26507; FLJ34646; GATF-C; GNAT; MGC15397; MGC15937 ;

**Immunogen:**

Human GLYATL1.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A