

GNAI3 Antibody



PACO09509

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GNAI3

Uniprot

P08754

Synonyms:

guanine nucleotide binding protein (G protein), alpha inhibiting activity polypeptide 3; GNAI3;87U6; FLJ26559 ;

Immunogen:

Human GNAI3.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A