

GNAQ Antibody



PACO09511

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GNAQ

Uniprot

P50148

Synonyms:

guanine nucleotide binding protein (G protein), q polypeptide; GNAQ; G-ALPHA-q; GAQ ;

Immunogen:

Human GNAQ.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A