

GNG2 Antibody



PACO09519

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GNG2

Uniprot

P59768

Synonyms:

guanine nucleotide binding protein (G protein), gamma 2;

Immunogen:

Human GNG2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20°C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A