

GNG4 Antibody



PACO09520

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GNG4

Uniprot

P50150

Synonyms:

guanine nucleotide binding protein (G protein), gamma 4; GNG4; DKFZp547K1018; FLJ23803; FLJ34187 ;

Immunogen:

Human GNG4.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A