

GPBP1 Antibody



PACO09543

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GPBP1

Uniprot

Q86WP2

Synonyms:

GC-rich promoter binding protein 1; GPBP1; DKFZp761C169; GPBP; MGC126339; VASCULIN ;

Immunogen:

Human GPBP1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A