

GPD1 Antibody



PACO09549

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GPD1

Uniprot

P21695

Synonyms:

glycerol-3-phosphate dehydrogenase 1 (soluble); GPD1; FLJ26652 ;

Immunogen:

Human GPD1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A