

GPR132 Antibody



PACO09558

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GPR132

Uniprot

Q9UNW8

Synonyms:

G protein-coupled receptor 132; GPR132; G2A; MGC99642 ;

Immunogen:

Human GPR132.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A