

# GPR137B Antibody



PACO09560

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

GPR137B

**Uniprot**

O60478

**Synonyms:**

G protein-coupled receptor 137B; GPR137B; TM7SF1 ;

**Immunogen:**

Human GPR137B.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A