

GPR161 Antibody



PACO09562

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:

Protein Background:

Gene ID:

GPR161

Uniprot

Q8N6U8

Synonyms:

G protein-coupled receptor 161; GPR161; FLJ33952; RE2 ;

Immunogen:

Human GPR161.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A