

GPR162 Antibody



PACO09563

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GPR162

Uniprot

Q16538

Synonyms:

G protein-coupled receptor 162; GPR162; A-2; FLJ78837; GRCA ;

Immunogen:

Human GPR162.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A