

GPR22 Antibody



PACO09567

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB

Recommended dilutions:**Protein Background:****Gene ID:**

GPR22

Uniprot

Q99680

Synonyms:

G protein-coupled receptor 22; GPR22; MGC129847 ;

Immunogen:

Human GPR22.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A