

GPS1 Antibody



PACO09571

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:**Protein Background:****Gene ID:**

GPS1

Uniprot

Q13098

Synonyms:

G protein pathway suppressor 1; GPS1; COPS1; CSN1; MGC71287 ;

Immunogen:

Human GPS1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A