

GSTCD Antibody



PACO09613

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GSTCD

Uniprot

Q8NEC7

Synonyms:

glutathione S-transferase, C-terminal domain containing; GSTCD; DKFZp686l10174; FLJ13273 ;

Immunogen:

Human GSTCD.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A