

GSTK1 Antibody



PACO09614

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

GSTK1

Uniprot

Q9Y2Q3

Synonyms:

glutathione S-transferase κ 1; GSTK1; GST; GST13; GST13-13; GSTK1-1 ;

Immunogen:

Human GSTK1.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A