

GTDC1 Antibody



PACO09626

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:

Protein Background:

Gene ID:

GTDC1

Uniprot

Q4AE62

Synonyms:

glycosyltransferase-like domain containing 1; GTDC1; mat-Xa ;

Immunogen:

Human GTDC1.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A