

# HAPLN4 Antibody



PACO09670

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse, Rat

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB, IHC

**Recommended dilutions:****Protein Background:****Gene ID:**

HAPLN4

**Uniprot**

Q86UW8

**Synonyms:**

hyaluronan and proteoglycan link protein 4; HAPLN4; BRAL2; KIAA1926 ; Brain link protein 2

**Immunogen:**

Human HAPLN4.

**Storage:**

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A