

HDHC2 Antibody



PACO09697

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC

Recommended dilutions:

Protein Background:

Gene ID:

HDHC2

Uniprot

Q7Z4H3

Synonyms:

HD domain containing 2; HDHC2; C6orf74; CGI-130; MGC87330; NS5ATP2; dJ167O5.2 ;

Immunogen:

Human HDHC2.

Storage:

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A