

# HDHD3 Antibody



PACO09699

---

## Product Information

**Size:**

50ul

**Reactivity:**

Human, Mouse

**Source:**

Rabbit

**Isotype:**

IgG

**Applications:**

ELISA, WB

**Recommended dilutions:**

**Protein Background:**

**Gene ID:**

HDHD3

**Uniprot**

Q9BSH5

**Synonyms:**

haloacid, dehalogenase-like hydrolase domain containing 3; HDHD3;2810435D12Rik; C9orf158; MGC12904 ;

**Immunogen:**

Human HDHD3.

**Storage:**

PBS with 0.1% Sodium Azide, 50% Glycerol, pH 7.3. -20&deg; C, Avoid freeze / thaw cycles.

## Product Images

---

N/A

N/A