

HEPACAM Antibody



PACO09707

Product Information

Size:

50ul

Reactivity:

Human, Mouse, Rat

Source:

Rabbit

Isotype:

IgG

Applications:

ELISA, WB, IHC, IF

Recommended dilutions:**Protein Background:****Gene ID:**

HEPACAM

Uniprot

Q14CZ8

Synonyms:

hepatocyte cell adhesion molecule; HEPACAM; FLJ25530; GlialCAM ;

Immunogen:

Human HEPACAM.

Storage:

PBS with 0.02% Sodium Azide, 50% Glycerol, pH 7.3. -20° C, Avoid freeze / thaw cycles.

Product Images

N/A

N/A